

OMRON Battery-operated Electronic Fever Thermometers

Applicable models: Refer to" EMC product reference table Electronic Fever Thermometers" Information for accompanying documents in the scope of EN60601-1-2:2015

Important information regarding Electro Magnetic Compatibility (EMC)
This device manufactured by OMRON HEALTHCARE Co., Ltd. conforms to EN60601-1-2:2015 Electro Magnetic Compatibility (EMC) standard. Nevertheless, special precautions need to be observed:
The use of accessories and cables other than those specified or provided by OMR

- The use of accessories and cables other than those specified or provided by OMRON could result in increased electromagnetic emission or decreased electromagnetic immunity of the device and result in improper operation.
- During measurement, the use of the device adjacent to or stacked with other device should be avoided because it could result in improper operation. In case such use is necessary, the device and other device should be observed to verify that they are operating normally.
- During measurement, portable RF communications device (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by OMRON. Otherwise, degradation of the performance of the device could result.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

Phenomenon	EMISSION Limits	Compliance			
Conducted and radiated RF EMISSIONS	CISPR 11	Group1, Class B			
Note: EMISSION tests for Harmonic distortion, Voltage fluctuations and flicker are not applicable.					

Table 1 - EMISSION Limits and Compliance



Table 2 - IMMUNITY TEST LEVELS

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS			
Electrostatic	IEC 61000-4-2	±8 kV contact			
discharge		±2 kV,±4 kV,±8 kV,±15 kV air			
Radiated RF	IEC 61000-4-3	10 V/m			
electromagnetic		80 MHz to 2.7 GHz			
fields		80 % AM at 1 kHz			
Proximity fields	IEC 61000-4-3				
from RF wireless communications		See table 3			
equipment					
Rated power	IEC 61000-4-8	30 A/m			
frequency		50 Hz and 60 Hz			
magnetic fields					
Note: IMMUNITY	lote: IMMUNITY tests for Input A.C. power port, Input D.C. power port, Patient coupling				
	ort and Signal input/output parts port are not applicable.				
For patient coupli	r patient coupling port is applicable to the electrostatic discharge.				

OMRON

Table 3 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications device

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM 土 5 kHz deviation 1 kHz sine	2	0.3	28
710			Pulse modulation 217 Hz	0.2	0.3	9
745	704 to 787	LTE Band 13, 17				
780						
810		GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720		GSM 1800; CDMA 1900;		2	0.3	28
1845		GSM 1900; Puls DECT;	Pulse modulation 217 Hz			
1970		LTE Band 1, 3, 4, 25; UMTS				
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n , RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5240 5500 5100 to 5800 5785	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						